
# Registration Matrices CT REF 55749-01, 55749-02, 55749-03

## General Information

**Registration Matrices CT** are designed for automatic patient registration in navigated spinal and cranial surgery. They are useful for intraoperative imaging with 3D Computed Tomography and Cone Beam CT.

**Registration Matrices CT** are for use with Brainlab navigation systems in combination with the Automatic Registration Software.

## Components

There are three different versions for different intended uses:

| No. | Component                              | Article Number |
|-----|--|----------------|
| ①   | Open Surgery (Spine)                   | 55749-01       |
| ②   | Small Incision (Spine)                 | 55749-02       |
| ③   | Minimally Invasive (Cranial and Spine) | 55749-03       |

## How to Use

**Registration Matrices CT** are designed for automatic patient registration. The detailed workflow is described in the **Automatic Registration Software User Guide**.

The expected lifetime of the device is 500 cycles of reprocessing and use.

Function can be ensured within the expected lifetime only. Continued use thereafter may increase the likelihood of failure or malfunction (inaccuracy).

Always attach new, unused marker spheres to the Registration Matrices prior to use.

Plan the OR setup prior to surgery. The camera must have an unobstructed view of the clean reflective marker spheres. The selected Matrix shall be positioned close to the region of interest and within the reconstructed 3D-volume of the scanner (for further information see the Automatic Registration Software User Guide).

Information on how to correctly attach marker spheres to the pins is available in the **Instrument User Guide**.

The matrices are highly accurate and sensitive medical devices. Handle them with care. If you drop or otherwise damage one, you must verify its calibration accuracy or return it immediately to Brainlab for testing, replacement or repair. Failure to do so may lead to serious patient injury.

Only use a **Brainlab Patient Scan Drape** to avoid the risk of contamination and to ensure that the marker spheres are visible to the camera.

Verify the registration before each surgical treatment to avoid patient injury.

## Reprocessing - General Information

**Must be reprocessed/sterilized before use. All components are delivered unsterile, unless specifically stated otherwise. Follow the reprocessing instructions prior to initial use and after each subsequent use.**

Due to the potential danger of residue, the application of rinsing aids is not recommended.

## Reprocessing Steps

| Steps |   |
|-------|---|
| 1.    | Pretreatment.   |
| 2.    | Automatic cleaning and disinfection.  |
| 3.    | Steam Sterilization (132 °C for at least five minutes): Store the Registration Matrices in the appropriate Sterilization Tray (52324) the right way around. Check the outlines on the bottom of the tray. |

Ensure that each Matrix fits correctly and doesn't protrude outside the tray. The lid must not touch any Matrix.

After cleaning and disinfecting, ensure that the Registration Matrices are clean and undamaged. Prior to use, verify that the devices are in good working condition. Do not use them if they are damaged or bent.

Repeated reprocessing, up to 500 cycles, has shown to have minimal effect on the device. End of life is normally determined by wear, damage and inaccuracy due to use. When a surgical device reaches the end of its functional life, clean it of all biomaterial/biohazards and safely discard the device in accordance with applicable laws and regulations.

Follow the service instructions in the Automatic Registration Software User Guide and perform the service checks as described on a regular basis to recognize potential damages.

## Further Information

For further information, contact Brainlab support or refer to the Brainlab user guides:

- **Automatic Registration Software User Guide**
- **Cleaning, Disinfection and Sterilization Guide**
- **Spine & Trauma Instrument User Guide**, including information on **Disposable Reflective Marker Spheres**.

This product may be covered by one or more patents or pending patent applications. For details, see:

<https://www.brainlab.com/patent/>

## Symbols

Caution

Warning

CE label shows that the Brainlab product complies with the essential requirements of European Medical Device Directive (MDD).

U.S. federal law restricts this device to sale by or on the order of a physician.

Reference / article number

Quantity of included items

Non-sterile


Consult instructions for use

Product manufactured with components that do not contain natural rubber latex.

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