Supporting Arm for Registration Matrix CT REF 53089

Fig./Abb. 1


Fig./Abb. 2


Fig./Abb. 3

Supporting Arm for Registration Matrix CT REF 53089

ENGLISH

General Information

Supporting Arm for Registration Matrix CT is designed for use with automatic registration software in navigated spinal and cranial surgery. The supporting arm can be directly attached to the OR table rail with the **Registration Matrix CT Cranial and Spine (Minimally Invasive)** (55749-03) positioned near the region of interest.

NOTE: The supporting arm is to be used unsterile and above the patient scan drape and is not autoclavable.

Components (Fig. 1)

Supporting Arm for Registration Matrix CT consists of the following components:

No.	Component	Article No.
①	OR Side Rail Adapter for Supporting Arm	53089-01
②	Supporting Arm for Registration Matrix CT	53089-02 & 53089-03

How to Use (Fig. 2)

The supporting arm is designed to position the **Registration Matrix CT**. The detailed workflow is described in the **Automatic Registration Software User Guide**.

Plan the OR setup prior to surgery. The camera must have an unobstructed view of the reflective marker spheres and the 3D volume must contain the blue body of the matrix in order to register automatically with the Registration Matrix CT.

NOTE: The supporting arm is a highly accurate and sensitive medical device. Handle it with care.

NOTE: If you drop or otherwise damage the product, you must verify its functionality or return it immediately to Brainlab for testing, replacement or repair. Failure to do so may lead to serious patient injury.

*NOTE: To avoid the risk of contamination and to ensure that the marker spheres are visible for the camera, use **Brainlab Patient Scan Drape** only.*

*NOTE: Follow the service instructions in the **Automatic Registration Software User Guide** and perform the service checks as described on a regular basis to recognize potential damage.*

Reprocessing (Fig. 3)

Step
1. Open the lever of the OR Side Rail Adapter for Supporting Arm and extend the telescope to its maximum position.
2. Remove any dirt and stains from all surfaces and gaps using a soft tissue moistened with surface disinfectants.

Do not use automatic cleaning and disinfection procedures.

Do not sterilize. Sterilization destroys the device and could damage sterilization equipment.

Compatible Disinfectants

Disinfectant Type	Example
Alcohol-based	Mikrozyd AF Liquid
Alkylamine-based	Incidin Plus 2%
Active oxygen based	Perform
Aldehyde/chloride based	Antiseptica Combi Surface Disinfection

Always closely follow the directions of the disinfectant manufacturer.

Using cleaning fluids, disinfection wipes or cleaning procedures other than those specified may damage equipment. It is recommended to only use disinfectants verified by Brainlab to avoid damaging the supporting arm.

Symbols

	Warning
	Caution
	The product complies with the essential requirements of the Medical Device Directive.
	U.S. federal law restricts this device to sale by or on the order of a physician
	Reference/article number

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