


# iPLAN NET WORKFLOW

Software Versions: iPlan Net 3.5




## STEP 1

- Open your browser.
- Enter the **iPlan Net** URL (e.g., https://iplan.net).
- In the **Login** tab, enter your user name and password and click the **Login** icon.
- You are now connected to the **iPlan Net** portal.



**NOTE:** To access additional documentation, press the **Help** button at the top of the **iPlan Net Connection Portal**.



## STEP 2A (IF APPLICABLE)

- If data transfer is not integrated in your **iPlan** version, select the relevant **PatXfer** application, and proceed through the software to convert patient data from DICOM to Brainlab format.
- To exit **PatXfer** correctly, press **Exit** followed by **OK** in the dialog which appears.

**NOTE:** As only one user may use **PatXfer** at a time, exit **PatXfer** correctly after successful data transfer.



## STEP 2B (IF APPLICABLE)


- Insert a USB flash drive into the USB port on the computer **before** opening **iPlan**.

**NOTE:** If you will export to a network location, you can skip this step and proceed to step 3.



## STEP 3


- Start **iPlan** by clicking on the relevant **iPlan** application button in the **Applications** tab.
- Perform relevant preoperative planning.



#### STEP 4 (USING INTERNET EXPLORER)

- Select your export location in the **iPlan** software and export your planned data.
- Wait for the data to transfer successfully.

**NOTE:** USB export is only possible with Internet Explorer. For HTML5-based connections using browsers other than Internet Explorer, you should export via the network to the end treatment device or storage location.



#### STEP 5

- Exit **iPlan** by selecting **Exit** in the **iPlan Navigator** (do not click X to close the window).
- Remove the USB flash drive from your computer (if applicable).
- Your planned data can now be used for patient treatment.


# SESSION SHARING

Software Versions: iPlan Net 3.5



## STEP 1


- Click **Session Sharing** in the Applications tab.



## STEP 2

- The **Session Sharing** dialog opens listing all active sessions.
- Select the required session and click **Share Session**.


**NOTE:** Advanced 3D sessions are not listed and cannot be shared.



## STEP 3

- To accept session sharing, the session owner must click **Yes** in the authorization notification dialog.
- The selected session is opened.

**NOTE:** Both the session owner and the guest user can now use the software's full functionality and see all steps performed by the other user.



## STEP 4


- To leave the shared session, click the **Disconnect** button at the bottom right of the screen.

**NOTE:** When you disconnect, the application does not close. You will simply leave the shared session.



# TABLET CONNECTIVITY

Software Versions: iPlan Net 3.5



## LOG-IN


- Open your HTML5 compatible browser.
- Navigate to the **iPlan Net** Server using the web browser on the tablet device.
- Enter the iPlan Net URL (e.g., <https://iplanet/help.aspx>).

**NOTE:** Wireless access to the server required.



## SELECT APPLICATION

- Start **iPlan** by clicking on the relevant **iPlan** application button in the **Applications** tab.
- Perform relevant preoperative planning.




## VIEWING

- It is recommended to use full screen capable browsers.

**NOTE:** Exact viewing functionality is browser-dependent and the chosen application may not always fit exactly to the size of the screen.

**NOTE:** Navigating back and forward should not be done during surgery planning.



## ACCESSING THE KEYBOARD

- Activate the keyboard by selecting the keyboard icon from the overlay menu.

**NOTE:** Exact functionality is browser-dependent.

**MANUFACTURER INFORMATION:**

Brainlab AG  
Olof-Palme-Str. 9, 81829 Munich, Germany

Europe, Africa, Asia, Australia: +49 89 99 15 68 44  
USA and Canada: +1 800 597 5911  
Japan: +81 3 3769 6900  
Latin America: +55 11 33 55 33 70  
France: +33-800-67-60-30

E-mail: [support@brainlab.com](mailto:support@brainlab.com)

**COPYRIGHT:**

This guide contains proprietary information protected by copyright. No part of this guide may be reproduced or translated without the express written permission of Brainlab.

Document Revision: 1.2  
Article Number: 60910-91EN

**LIABILITY:**

This guide is subject to change without notice and does not represent a commitment on the part of Brainlab.

For further information, please refer to the "Limitations of Liability" section in the Brainlab Standard Terms and Conditions of Sale.

