


CT Protocol – HIP

How to scan hip data for BrainLAB navigation:

Patient orientation	<ul style="list-style-type: none"> • Strict supine position, • Legs parallel with no internal or external rotation • Orientation “Feet first” recommended
<p>Slice Thickness / Scan Range</p>  <p><i>Range 1: Pelvis</i></p> <p><i>Range 2: Femur</i></p> <p><i>Sample Scout View</i></p> <p><i>x, y remain constant</i></p>	<ul style="list-style-type: none"> • Include the whole pelvis and both upper femurs in the scan. • The amount of spine seen in the field of view should be kept to a minimum. • Include as little tissue as possible; only bone structure should be visible. • The width and height (x,y) of both FOVs must be equal, so that the center points have the same x, y coordinates. • Range 1: Entire Pelvis: Start scan at the cranial end of the iliac crest and scan to 10 cm below the lesser trochanter . The pelvis should be centered in the field of view. <i>Slice Thickness: 2-3mm</i> <p>Recommendation: helical slices at 5mm reconstructed to 2mm. Not more than 250 slices in the entire dataset</p> <p>Reconstructions are allowed!</p> <ul style="list-style-type: none"> • Range 2: Femoral Condyles: ca. 50mm of the left and right distal femurs <i>Slice Thickness: maximum 5mm</i>
Image / Pixel Size	<ul style="list-style-type: none"> • Pixel size must be the same for both Field of Views!
Scan Properties	<ul style="list-style-type: none"> • Scan with soft tissue windowing
Scan Technique	<ul style="list-style-type: none"> • Sequential scans: continuous or overlapping slices. No gap! • Helical Scans: pitch (table:scan ratio) =1:1 recommended. <i>Reconstructed images are allowed.</i> • Slice thickness can be changed during the scan.
Gantry tilt	<ul style="list-style-type: none"> • No gantry tilt allowed!
Table height	<ul style="list-style-type: none"> • Must remain the same during the scan.
Matrix Size	<ul style="list-style-type: none"> • Necessary: 512x512
Scan direction	<ul style="list-style-type: none"> • cranial to caudal • caudal to cranial
Storing	<ul style="list-style-type: none"> • Store both scanning sets as one patient file.

* Based on PatXfer 4.21 or higher and VV² Hip Navigation Software

* If you need additional information please contact your local BrainLAB Support Engineer